Observation

Microscopie

Observation

Microscopy

Réf :
571 248
571 249

Français – p 1

English – p 9

Version : 0104

Microscopes OPTIKA série B-190
Microscopes OPTIKA B-190 series

Sommaire

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0825 563 563 *

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Contents

1	Description	10
2	Introduction	12
3	Unpacking and assembly	12
4	Using the microscope	12
5	Maintenance	13
6	Technical specifications	14
7	Recovery and recycling	15
8	After-Sales Service	15

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
Power supply
input jack

2 Introduction

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

Jeulin reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Jeulin declines any responsibility deriving from instrument uses that do not comply with this manual.

3 Unpacking and assembly

3.1 Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

3.2 Assembly

Place the observation head onto the top of the arm and tighten the lock-screw. Insert the eyepieces into the eye tubes.

3.3 Plug-in

Connect the provided 7,5Vdc power supply to the power supply input jack on the rear of the microscope.

4 Using the microscope

4.1 Adjust the observation head

Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw.

4.2 Place the specimen on the stage

Lock the specimen slide on the mechanical stage using the slide clamp. Ensure that the specimen is centred over the stage opening by adjusting the coaxial knobs of the stage.

4.3 Lamp settings

The microscope is fitted with power white LED illuminator. Before turning on the light, read the section 5.3 about electrical safety precautions. Insert the plug of the 7,5Vdc power supply into the jack on the rear and turn on the ON/OFF switch. Turn the brightness adjustment knob to a brightness suitable for observation.

4.4 Adjust interpupillary distance

Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen.

4.5 Focus and diopter adjustment

Rotate the coarse focusing knob to bring the slide into focus using an objective with a low magnification. Adjust the fine focusing knob to get the image sharp and clear while observing with your right eye, then turn the dioptic adjustment ring on the left eyepiece until a sharp and clear image can be observed also with the other eye. Turn the tension-adjust-knob to get a suitable tension for the focus system. When the image is in focus, turn the nosepiece to choose the objective you need.

4.6 Condenser adjustment

Hold the notched ring at the base of the condenser and rotate it up or down until you reach a clear and uniform sample illumination.

4.7 Set the numerical aperture

Adjust the aperture of the iris diaphragm to set the numerical aperture of the illuminator, thus controlling image contrast and resolution.

4.8 Additional filters

The green glass filter can be inserted in the flip-out filter holder underneath the condenser.

5 Maintenance

5.1 Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85% (non condensing). Use a dehumidifier if needed.

5.2 To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

5.3 Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

5.4 Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.

Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.

- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.


5.5 Instructions to return your product to the supplier

If you need to send the microscope to Jeulin for maintenance, please use the original packaging.

6 Technical specifications

HEAD:	B-191:	monocular, 360° rotating, 30° inclined
	B-192:	binocular, rotating, inclined
EYEPIECES:		WF 10X/ 18mm
NOSEPIECE:		Quadruple, reversed
OBJECTIVES:		Achromatic 4x, 10x, 40x, 100x (oil immersion)
WORKING STAGE:		Double layer with mechanical sliding stage, 125x115mm, moving range 50x30mm
CONDENSER:		1.25N.A. Abbe type, with iris diaphragm
FOCUSING SYSTEM:		Coaxial coarse and fine focusing system, with tension adjustment mechanism
ILLUMINATION:		Extra-power White LED, non-rechargeable, with brightness control.
MAINS POWER:		Input Voltage of the power supply: 100-240Vac / 50-60Hz

7 Recovery and recycling

The appliance reports the symbol:

This symbol means that the appliance can be a precious source of raw materials. Therefore, it must not be disposed of as waste, but separately collected for the recycling and the recovery of the materials it contained in it. Such materials, if improperly dispersed into the environment, can be harmful to the environment and to human health.

The producer of the equipment, Jeulin, recovers, re-uses and recycles the raw materials contained in the equipment. Such recovery, however, needs your help.

When, at the end of its operating life, you shall decide to dispose of the apparatus, do not try to open it, nor to use parts of it in ways other than reported in this Manual, but bring it back to the Distributor. The Distributor shall collect the apparatus free of charge.

The recovery of the raw materials shall then take place in accordance with the European Directive 2002 / 96/ EC and all other relevant Directives. Never disassemble, nor dispose of as waste, apparatus reporting the "crossed bin" symbol indicated above.

8 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries